

Nirav 2 Section Treatment Table is meticulously designed with a short head section for precise manipulations. Easily adjustable, it allows optimal patient positioning and offers exceptional versatility for effective treatment of the neck and shoulder areas.



## KEY FEATURES



## DIFFERENCE BETWEEN PRO AND BASIC MODEL

### PRO MODEL

- **Walk Free Switch**

The Pro Model comes with a hands-free foot control system for height adjustments. It enhances hygiene, efficiency, ergonomics and can be operated from both sides. IEC 60601-1 safety requirement and IEC 60601-1-2 EMC standard complaint.

- **Easy Mobility**

The Pro Models have an easy mobility feature, which allows effortless repositioning of the treatment table within clinical environments. With a press of a lever, the table is lifted onto four wheels for smooth movement. When mobility is no longer required, the table can be securely fixed back on the ground for stability during procedures.

### BASIC MODEL

- **Foot Switch**

Basic Model only comes with a wired foot controller for height adjustment by pressing the + or – symbols on the foot switch, as shown in the image below.



## UPHOLSTERY

### Phthalate-Free Upholstery

Phthalate-Free Upholstery offers a safer, hospital-grade surface that resists stains, disinfectants, and abrasion, ensuring hygienic, durable, and eco-friendly performance across demanding clinical environments.



## FEATURES

|                  |  |
|------------------|--|
| Backrest         | Manually Adjustable (Gas Spring)                               |
| Nose Slot        | Helps the patient breathe freely in face down (prone position) |
| Height           | Electronically Adjustable                                      |
| Walk Free Switch | For Hands Free Height Adjustment                               |
| Easy Mobility    | 75 mm Castor Wheels for smooth movement                        |

## TECHNICAL SPECIFICATIONS

|                      |                   |
|----------------------|-------------------|
| Length               | 1980 mm           |
| Width                | 750 mm            |
| Minimum Height       | 530 mm            |
| Maximum Height       | 890 mm            |
| Backrest Inclination | 32° (when raised) |

## ELECTRICAL PARAMETERS

|                               |                         |
|-------------------------------|-------------------------|
| Power in (Mains Power Supply) | 100-240 V, 50-60 HZ, 2A |
| IPX Rating                    | IPX4                    |
| Operating Temperature         | 5-40°C                  |
| Shock Protection:             | Class I                 |

## BASE & FRAME

|                               |                          |
|-------------------------------|--------------------------|
| Base                          | Powder Coated Mild Steel |
| Safe Working Load             | 200 kgs                  |
| Patient Load Bearing capacity | 150 kgs                  |

## SAFETY STANDARDS

|   |             |
|---|-------------|
| General safety for medical electrical equipment | IEC 60601-1 |
| EMI/EMC   | IEC 60601-2 |

## UPHOLSTERY

|                  |                                    |
|------------------|------------------------------------|
| Type             | Phthalate-Free Upholstery          |
| Thickness        | 1.01 mm                            |
| Tensile Strength | ISO 1421, 804 N/ 5 cm              |
| Colour Options   | Cream, Beige, Brown, Mouse Grey    |
| Weight           | 732 g / sqmt                       |
| Abrasion         | ISO 5470-2 = 2,00,000 cycles / min |

## WARRANTY

3 years on motors

