

The Esthetica Basic is a stationary phlebotomy chair, built for the daily demands of clinics, diagnostic laboratories, blood banks, and hospital phlebotomy rooms. The patient's supporting arm rests on a padded flip armrest and elbow rest, both adjustable to hold the arm comfortably in place, while the dedicated blood drawing arm on the opposite side positions the working arm at the ideal height for venous access. Built on a heavy-duty mild steel sheet metal and tube frame and finished in durable, easy-clean phthalate-free upholstery, it delivers the stability and hygiene that high-frequency clinical use requires — in a fixed-height format that keeps things simple and dependable.



# KEY FEATURES

## Height Adjustable Blood Drawing Arm

On the draw side, the dedicated blood drawing arm positions the patient's arm at the ideal height for venous access. Height adjustable to suit patients of different builds, it helps achieve comfortable arm positioning while providing the phlebotomist with a stable working surface for blood draws and other phlebotomy procedures.

## Padded Flip Armrest with Height & Depth Adjustment

On the support side, a padded flip armrest holds the patient's non-draw arm in a comfortable, secure position. Adjustable in both height and depth to suit each patient, it flips clear of the seat when the session is done for easy, unassisted exit.

## Height Adjustable Padded Elbow Rest

Alongside the flip armrest, a height adjustable padded elbow rest gives the patient's supporting arm a stable resting point throughout the procedure — keeping them relaxed and still while the draw is carried out.

## Padded Seat & Backrest

The generously padded seat and backrest are designed to keep patients comfortable during blood draws, transfusions, and other clinical procedures. High-density cushioning helps reduce pressure points and supports correct seated posture, making the chair suitable for both short and extended treatments.

## Fixed Footrest


A dedicated fixed footrest provides comfortable lower-leg support while seated, helping patients maintain a stable, relaxed position throughout the procedure. Its integrated design adds comfort without introducing moving parts, contributing to the chair's long-term durability and reliability.

## Heavy-Duty Sheet Metal, Tubular Metal & Stainless Steel Frame

The Esthetica Basic is built on a robust frame combining sheet metal, tubular metal, and stainless steel components for strength, stability, and long service life. The fixed base structure utilizes heavy-duty sheet metal and tubular metal construction, while stainless steel pipes are incorporated in key functional elements. Finished with a durable powder-coated coating, the frame is designed to withstand the demands of high-frequency clinical use.

## Phthalate-Free Upholstery

The Esthetica Basic is upholstered in phthalate-free upholstery, designed for professional clinical environments. It resists stains, oils, and disinfectants, withstands frequent cleaning, and carries antibacterial and fire retardant properties as standard..



Cream

Beige

Brown

Mouse Grey