ESTHETICA
Medical Furniture

DIALYSIS CHAIR
4 MOTOR MODEL

The Esthetica Dialysis Chair is designed to provide dependable comfort, safety, and support for patients during long haemodialysis sessions while ensuring smooth and reliable operation in clinical environments. With ergonomic positioning, robust construction, and user-friendly controls, it supports efficient clinical workflows while maintaining optimal patient comfort throughout treatment.



KEY FEATURES



Electronic Footrest Adjustment

Allowing nursing staff to position patients easily with the push of a button. This eliminates the need for manual force required in conventional chairs, reducing physical strain on caregivers and enabling faster, smoother adjustments during treatment.

Two Pre-Set Positions & Quick Release Function

End of Treatment and Trendelenburg and a quick-release backrest enable medical staff to rapidly adjust the chair according to clinical needs, improving workflow efficiency and ensuring patient safety during procedures or emergencies.

Lifting Column

A 3-segment lifting column provides greater height adjustment, improved stability, and smoother, quieter motion than a linear actuator. Its telescopic design saves space while offering higher lifting capacity and better static and bending load performance.

Phthalate-Free Upholstery

Phthalate-Free Upholstery is designed for hospital environments, offering a safer surface for patients and healthcare workers. It resists stains and disinfectants, withstands frequent cleaning, and ensures durable, fire- and abrasion-resistant performance.


End Of Treatment



CPR


Cream Beige Brown Mouse Grey

FEATURES

Armrests	Adjustable Armrest with Back & Forth and Up - Down Movements
Backrest	Electronically Adjustable
Legrest	Electronically Adjustable
Foot Rest	Electronically Adjustable
Height	Electronically Adjustable
Trendelenburg	Integrated with Legrest
Castors	100 mm high grade synthetic, single wheel castos, four with brakes
Memory Positions	End of Treatment; CPR

TECHNICAL SPECIFICATIONS

In Flat Bed Position:

Length	2010 mm / 74.4 in
Width without side armrests	570 mm / 22.4 in
Width with side armrests	880 mm / 34.6 in
Minimum seat height	550 mm / 21.6 in
Maximum seat height	700 mm / 27.5 in
Minimum height	1180 mm / 46.4 in
Maximum height	1330 mm / 52.3 in

ELECTRICAL PARAMETERS

Power in (Mains Power Supply)	100–240V, 50–60Hz 1.5A
IPX Rating	IPX6
Operating Temperature	10–40°C
Shock Protection:	Class I

ADD-ON ACCESSORIES

Stainless steel IV stand
ABS food tray
IPad/Mobile Stand — Adjustable height and angle, positionable on either side
Dual power socket panel for charging laptops, tablets, and phones

BASE & FRAME

Base	Powder-coated mild steel frame
Safe Working Load	150 kgs
Patient Load Bearing capacity	135 kgs

CERTIFICATIONS & STANDARDS

Quality management	ISO 13485:2016 certified manufacturer (MTIC)
General safety for medical electrical equipment	IEC 60601-1
EMI/EMC	IEC 60601-1-2


UPHOLSTERY

Type	Veridura Phthalate-Free Upholstery
Thickness	1 mm
Tensile Strength	ISO 1421, 869.22 N/ 5 cm
Colour Options	Cream, Beige, Brown, Mouse Grey
Weight	726.7 g / sqmt
Abrasion	ISO 5470-2 = 2,00,000 cycles / min

Features : Phthalate-Free; Fire Retardant; Oil Resistant; Alcohol Resistant; Ultra Soft; Stain Resistant; Anti Bacterial; Eco-Friendly

WARRANTY


3 years on motors





HANDSET

TECHNICAL DRAWINGS


FORNT VIEW








IN FLAT BED POSITION



MINIMUM HEIGHT



MAXIMUM HEIGHT

ALL DIMENSIONS ARE IN MM